


From consultation to recovery and restoration,
quickly, efficiently and confidently.



SIG

Scan Implant Guide



Smart Planning, Safe Placement, Simple Solution.

Call to Order 855.373.1614

5 Simple Steps

Our simple, five-step process allows the doctor to take the patient from initial consultation to recovery and restoration, quickly, efficiently and confidently.

1. Scan the Patient

Take a Cone Beam CT scan of the patient (or, alternately, scan the patient wearing a radiographic template and scan the template separately).




3. Plan the Surgery

Upload the scans to our Scan Implant Guide Planning Software and plan the implant surgery. You can either acquire the software and plan the surgery yourself, or have one of our experienced clinicians assist you.



2. Scan the Model

Make a stone model of the patient's arch and send it to us for 3D scanning (an intra-oral surface scan of the patient is also acceptable).



4. Print the Guide

Finalize your order, approve the surgical plan and send it to us to fabricate the surgical guide with our state-of-the-art 3D printer.



5. Place the Implants

Perform the implant surgery with confidence, using a precision-crafted Scan Implant Surgical Guide.



SIG
Scan Implant Guide

2999 Overland Ave. #116
Los Angeles, CA 90064
Tel: 855.373.1614
Fax: 310.730.6241
www.si-guide.com

Learn how we can help you simplify implant surgery.